


- Concept

The concept

When small ice cubes are added, submerged and stirred within a calorimeter with water, the heat gained by the ice cubes to melt and obtain the equilibrium temperature will be equal to the heat lost by calorimeter, stirrer and water.


- Initial Setup



- Measurements you need to take in the correct order.

1. Measure the mass of the calorimeter with the stirrer. $M_1 - (M_c)$
2. Fill the calorimeter with water and then measure the mass. $M_2 - (M_c + M_w)$
3. Initial temperature of water θ_1
4. Get the lowest temperature after adding ice θ_2
5. Final mass of calorimeter + stirrer + ice M_3

(Here $\theta_1 > \theta_2$



- A - Temperature of calorimeter + stirrer + water decreasing when ice cubes are added.
- B - 0°C ice turning into 0°C water.
- C - 0°C water getting heated and coming into an equilibrium with water in the calorimeter.

After taking these readings, we can find the specific heat capacity of the ice using the equation $H = ML$ for the phase change of ice to water and we use $H = ms\theta$ for the other heat exchanges

So we assume that the heat gained from the ice cubes are equal to the heat released by the water and calorimeter

$M_s = \text{Mass of ice added } (M_3 - M_2) = \text{Mass of water added}$
 (As the same ice is melt in to water)

$M_w = \text{Mass of initial water in calorimeter } (M_2 - M_1)$

$M_c = \text{Mass of the calorimeter + Stirrer } (M_1)$

$L = \text{Spefic latemt heat of fussion of water}$

$S_w = \text{Spefic heat capacity of water } (4200)$

$S_c = \text{Spefic heat capacity of calorimeter } (4000)$

$$H = mS\delta\theta = ML$$

$$M_s L + M_s S_w (\theta_2 - 0) = M_c S_c (\theta_1 - \theta_2) + M_w S_w (\theta_1 - \theta_2)$$

$$\therefore L = \frac{(M_c S_c + M_w S_w)(\theta_1 - \theta_2) - M_s S_w (\theta_2 - 0)}{M_s}$$

Important point

- Why should we use a net stirrer to stire this instead of a normal one?

As the density of ice is less than water ice floats on water. Therefore to keep them in the bottom we need a net stirrer



- Why should we wipe the ice cubes with a blotting paper before adding it to the calorimeter?

To make sure its dry and no water is added with the ice cubes

- Which shape of ice is best for this experiment?

small-medium cubes

Large onces can't be used as we can't mantain the temperature the calorimeter gains and it might have a temperature difference between the surface and the inner temperature (inner

